

Nucleosome R-MetStat Panel

Cat. No. NUC-068 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	A 96-well plate containing 15 unique mononucleosomes assembled from recombinant human histones expressed in E. coli.
Source	E.coli
Description	A 96-well plate containing 15 unique mononucleosomes assembled from recombinant human histones expressed in E. coli (two each of histones H2A, H2B, H3 and H4; accession numbers: H2A-P04908, H2B-O60814, H3.1-P68431/H3.2-Q71DI3, H4-P62805). Single histone post-translational modifications (PTMs) are created using a proprietary synthetic method. The nucleosomes are wrapped by 147 base pairs DNA containing a central 601-positioning sequence (Lowary and Widom 1998 J Mol Biol. (1):19-42) with a 5' biotin-TEG group.
Form	Purified recombinant mononucleosomes individually stored in a 96 well plate (1.5 µg nucleosome at 1.5 µM per well) in 5 µL of nucleosome storage buffer (10 mM Tris-HCl pH 7.5, 1 mM EDTA, 25 mM NaCl, 2 mM DTT, & 20% glycerol).
Molecular Mass	~200,000 Da
Applications	The access to epigenetic diversity in a physiological nucleosome context enables broad end-user applications, including nucleosome binding studies (e.g. chromatin reader binding preferences), enzyme screening assays (e.g. identification of preferred substrates), and antibody specificity testing (e.g. for applications that require antibody specificity in a nucleosome, rather than histone peptide context). The biotin group on

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

the DNA facilitates pull-down experiments. Mononucleosomes do not contain free DNA which could alter assayed activities.

Storage

Stable for six months at -80°C from date of receipt.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA