

Human Melanoma LOX-IMVI Whole Cell Lysate

Cat. No. LOX-IMVI-060WCY **Lot. No.** (See product label)

SPECIFICATION

Product Overview	LOX-IMVI (Human Melanoma) Whole Cell Lysate
Species	Human
Source	Skin
Description	LOX-IMVI human melanoma cell line is an excellent model for probing mechanisms of metastasis and for evaluation of chemotherapies.
Applications	IP, WB
Stability	The lysate is stable for at least one year from date of shipping when stored at -80 centigrade. After incubation with SDS-PAGE sample buffer, SDS-PAGE samples need to be stored in -20 centigrade for long term storage.
Storage	Store the lysate at -80 centigrade and avoid repeated freeze / thaw cycles.
Concentration	1 vial of 100 µg whole cell lysate in RIPA buffer. 1 vial of 250 µL 2 x SDS Sample Buffer (4% SDS, 125 mM Tris-HCl pH 6.8, 10 % Glycerol, 0.002 % Bromphenol blue, 100 mM DTT)
Storage Buffer	Modified RIPA buffer (25 mM Tris-HCl pH 7.6, 150 mM NaCl, 1% NP-40, 1 mM EDTA, 1 x Proteinase inhibitor cocktail mix (Sigma) and 1 mM Na ₃ VO ₄)
Shipping	Shipped on dry ice.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA