

Recombinant *Bacillus subtilis* subsp. Natto aprN protein, His-tagged

Cat. No. aprN-183N **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Recombinant <i>Bacillus subtilis</i> subsp. Natto aprN protein(P35835)(107-381aa), fused with N-terminal His tag, was expressed in <i>E. coli</i> .
Species	<i>Bacillus subtilis</i>
Source	<i>E. coli</i>
ProteinLength	107-381aa
Form	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Molecular Mass	31.8 kDa
AA Sequence	AQSVPYGISQIKAPALHSQGYTGSNVKVAVIDSGIDSSHPDLNVRGGASFVPSETNP YQDGSSHGTHVAGTIAALNNSIGVLGVAPSASLYAVKVLDTGSGQYSWIINGIEWAI SNNMDVINMSLGGPTGSTALKTVVDKAVSSGIVVAAAAGNEGSSGSTSTVGYPKY PSTIAVGAVNSSNQRFSSVGSELDVMAPGVSIQSTLPGGTYGAYNGTSMATPHV AGAAALILSKHPTWTNAQVRDRLESTATYLGNSFYFGKGLINVQAAAQ
Purity	Greater than 85% as determined by SDS-PAGE.
Storage	Store at -20°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

freeze-thaw cycles.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA