

Native Soybean Glycine max Lectin (SBA), Colloidal Gold-Conjugated

Cat. No. SBA-110S **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Native Soybean Glycine max Lectin (SBA) - Colloidal Gold Blood Group Specificity: Non-specific Preferred Sugar Specificity: α-N-Acetylgalactosamine > β-N-Acetylgalactosamine Inhibiting or Eluting Sugar: N-Acetylgalactosamine Divalent Ions: Ca⁺⁺, Mn⁺⁺ Excitation: 550 nm Emission: 570 nm Mitogenic Activity: Yes Hazardous Shipping: Non-hazardous
Species	Soybean
Source	Soybean
Description	Colloidal gold can be used to visualize localization of specific saccharide moieties and examine the distribution of plasma membrane glycoproteins based on the binding pattern of Glycine max (Soybean) Lectin (SBA). Colloidal gold labeled lectins can be viewed via electron microscopy.
Form	Liquid
Purity	High Purity Grade
Stability	Shelf Life: 1 year

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Storage	At 2-8 centigrade.
----------------	--------------------

Concentration	1 mg/mL
----------------------	---------

GENE INFORMATION

Official Symbol	SBA
------------------------	-----

Synonyms	Glycine max Lectin; SBA
-----------------	-------------------------

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA