

Native Glycine max Lectin (SBA), Separopore® 4B OnePASS™ Column-Conjugated

Cat. No. SBA-628S Lot. No. (See product label)

SPECIFICATION

Product Overview	Native Soybean Glycine max Lectin (SBA) - Separopore® 4B OnePASS™ Column Hazardous Shipping: Non-hazardous
Species	Soybean
Source	Soybean
Description	Glycine max agglutinin (SBA) is isolated and purified by affinity chromatography. It consists of four subunits of identical size. SBA has specificity for blood groups A1, A2, and B, and interacts better with neuramidase-treated cells than untreated cells. This lectin has selective affinity for lymphocytes and human CD34+ hematopoietic stem cells. Immobilized conjugates of SBA are important tools for removing T-cells in bone marrow transplants. Glycine max Lectin (SBA) OnePASS™ Spin Columns are standard SBA matrix prepacked in a column which can be used via gravity flow or as a spin column to rapidly separate bound and passthrough materials.
Storage	At 2-8 centigrade.

GENE INFORMATION

Official Symbol	SBA
Synonyms	Glycine max Lectin; SBA

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA