

Native Canine Collagen Type I and III Protein

Cat. No. Collagen Type I & III-05C **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Purified type I & III canine collagen.
Species	Canine
Source	Canine
Description	Type I collagen is the most abundant collagen and is found in connective tissues including tendon, ligament, dermis and blood vessel. It is the major component and the primary determinant of tensile strength of the extracellular matrix (ECM). It is widely used as a thin layer on tissue-culture surfaces to enhance the attachment and proliferation of a variety of cells including endothelial cells, fibroblasts, hepatocytes, epithelial cells etc. In addition, collagen I can self-assemble into a 3-D superamolecular gel in vitro, making it an ideal biological scaffold to promote more in vivo-like cellular morphology and function. Type III collagen is the second most abundant collagen in tissues and is found most commonly in tissues exhibiting elastic properties such as skin, lungs, intestinal walls and walls of blood vessels. It is a homotrimer comprised of three alpha-1 chains and resembles other fibrillar collagens in structure and function. It is synthesized as procollagen, similar to collagen I, but the N-terminal propeptide remains attached in the mature fibrillar type III form.
Form	Lyophilized, salt free, 10 mg/vial
Purity	Canine collagen type I and III - 98%

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Canine Collagen (other types) - <1%
Non-collagen proteins - below 0.5%

Storage

Collagen dissolved in acetic acid is stable at 4 centigrade for 1 month. Lyophilized collagen long term storage (2 years) at -20 centigrade or lower.

Reconstitution

Use 0.5 M acetic acid, pH 2.5. Dissolved collagen retains immunologic properties of native collagen. Structure of native collagen confirmed by ability to form microfibrils.

GENE INFORMATION**Official Symbol**

Collagen Type I and III

Synonyms

COL1A1; COL1A2; COL3A1; Collagen type III, alpha 1; osteogenesis imperfecta;
Collagen Type I and III

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA