

Recombinant Rat tapeworm Dihydrolipoamide dehydrogenase protein, His&Myc-tagged

Cat. No. Dihydrolipoamide dehydrogenase-6009R **Lot. No.** (See product label)

SPECIFICATION


Product Overview	Recombinant Rat tapeworm Dihydrolipoamide dehydrogenase protein(P80647)(1-53aa), fused to N-terminal His tag and C-terminal Myc tag, was expressed in E. coli.
Species	Rat tapeworm
Source	E.coli
ProteinLength	1-53aa
Form	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Molecular Mass	10.3 kDa
Purity	Greater than 90% as determined by SDS-PAGE.
Storage	Store at -20°C/-81°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.
Reconstitution	Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL.We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 51%.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

SDS-PAGE



Tel: 1-631-559-9269 1-516-512-3133

Email: info@creative-biomart.com Fax: 1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA