

# Recombinant Full Length *Oenothera Argillicola* Photosystem II Cp47 Chlorophyll Apoprotein(Psbb) Protein, His-Tagged

**Cat. No.** RFL30494OF    **Lot. No.** (See product label)

## SPECIFICATION

|                         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b> | Recombinant Full Length <i>Oenothera argillicola</i> Photosystem II CP47 chlorophyll apoprotein(psbB) Protein (P19819) (1-508aa), fused to N-terminal His tag, was expressed in <i>E. coli</i> .                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Species</b>          | <i>Oenothera argillicola</i> (Appalachian evening primrose)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Source</b>           | <i>E. coli</i>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>ProteinLength</b>    | Full Length (1-508)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Form</b>             | Lyophilized powder                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>AA Sequence</b>      | MGLPWYRVHTVVLNDPGRLLAVHIMHTALVAGWAGSMALYELAVFDPSPDVLDP<br>WRQGM FVIPFMTRLGITNSWGGWSITGGTVTNPGIWSYEGVAGSHILFSGLCFLAAI<br>WHWVYWDL AIFSDERTGKPSLDLPKIFGIHLFLSGLACFGFGAFHVTGLYGPGIWVS<br>DPYGLTGEVQP VNPAWGVGEFDPFVPGGIASHHIAAGTLGILAGLFHLSVRPPQRL<br>YKGLRMGNIETVLSS SIAAVFFAAAFVAGTMWYGSATTPIELFGPTRYQWDQGYFQ<br>QEIYRRVGAGLAKNQSLSE AWSKIPEKLAFYDYIGNNPAKGGLFRAGSMDSGDZIA<br>VGWLGHPIFRDKEGRELFVRRMP TFFETFPVVLVDGDGIVRADVPFRRAESKYSVE<br>QVGVTVEFYGGELNGVSYSDPATVKKY ARRAQLGEIFELDRATLKSDGVFRSSPRG<br>WFTFGHASFALLFFFGHIWHGARTLFRDVFA GIDPDLDTQVEFGAFQKLGDPTTRR<br>QAV |
| <b>Purity</b>           | Greater than 90% as determined by SDS-PAGE.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |

 Tel: 1-631-559-9269    1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)     Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

|                         |                                                                                                                                                                                                                                                                                                                                                                                                             |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Applications</b>     | SDS-PAGE                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Notes</b>            | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Storage</b>          | Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.                                                                                                                                                                                                                                                                                              |
| <b>Storage Buffer</b>   | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Reconstitution</b>   | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>GENE INFORMATION</b> |                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Gene Name</b>        | psbB                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Synonyms</b>         | psbB; Photosystem II CP47 reaction center protein; PSII 47 kDa protein; Protein CP-47                                                                                                                                                                                                                                                                                                                       |
| <b>UniProt ID</b>       | <a href="#">P19819</a>                                                                                                                                                                                                                                                                                                                                                                                      |