

## Recombinant Crimean-Congo hemorrhagic fever virus GP protein, MBP&His-tagged

Cat. No. GP-1197C    Lot. No. (See product label)


### SPECIFICATION

|                         |                                                                                                                                                                                                                                    |
|-------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b> | Recombinant Crimean-Congo hemorrhagic fever virus GP protein(Q8JSZ3)(248-519aa), fused to N-terminal MBP tag and C-terminal His tag, was expressed in E. coli.                                                                     |
| <b>Species</b>          | Crimean-Congo hemorrhagic fever virus                                                                                                                                                                                              |
| <b>Source</b>           | E.coli                                                                                                                                                                                                                             |
| <b>ProteinLength</b>    | 248-519aa                                                                                                                                                                                                                          |
| <b>Form</b>             | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0. |
| <b>Molecular Mass</b>   | 74.5 kDa                                                                                                                                                                                                                           |
| <b>Purity</b>           | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                        |
| <b>Storage</b>          | Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.                                                                                                                     |
| <b>Reconstitution</b>   | It is recommended that sterile water be added to the vial to prepare a stock solution of 0.2 ug/ul. Centrifuge the vial at 4°C before opening to recover the entire contents.                                                      |

 Tel: 1-631-559-9269    1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)     Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

**SDS-PAGE**

Tel: 1-631-559-9269 1-516-512-3133

Email: [info@creative-biomart.com](mailto:info@creative-biomart.com) Fax: 1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA