

Recombinant Full Length *Korarchaeum Cryptofilum* Digeranylgeranylglyceryl Phosphate Synthase (Kcr_1480) Protein, His-Tagged

Cat. No. RFL26689KF **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Recombinant Full Length <i>Korarchaeum cryptofilum</i> Digeranylgeranylglyceryl phosphate synthase (Kcr_1480) Protein (B1L6Z7) (1-281aa), fused to N-terminal His tag, was expressed in <i>E. coli</i> .
Species	<i>Korarchaeum cryptofilum</i>
Source	<i>E. coli</i>
ProteinLength	Full Length (1-281)
Form	Lyophilized powder
AA Sequence	MGKLGAYVEISRPKNALMSILGTLTGWVNSTSVYDGRLLACLIPPLVLMAGNAINDY YD AEIDAINKPYPPIPSGRISKREALNIYIALSLFGIALSIFLGFIEFLIVTAFSLSWYAYA RWLKRTGVPGNALVSLGVAFTLIFGSLAAGNLTKNVIIFSSVAFTSNLIREFVKAVEDL P GDRAHGVRTIAVRIGVKRTGILVFLSLATVVLTLIPVIFRLTGIIYLSLSVIISLPILM L ASAICLKKGKLEERARETSSLIKVSMFLGLLGMLLDPFRVV
Purity	Greater than 90% as determined by SDS-PAGE.
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Storage	Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.
Storage Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
GENE INFORMATION	
Gene Name	Kcr_1480
Synonyms	Kcr_1480; Digeranylgeranylglyceryl phosphate synthase; DGGGP synthase; DGGGPS; (S-2,3-di-O-geranylgeranylglyceryl phosphate synthase; Geranylgeranylglycerol-phosphate geranylgeranyltransferase
UniProt ID	B1L6Z7

 Tel: 1-631-559-9269 1-516-512-3133 Email: info@creative-biomart.com  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA