

Recombinant Full Length *Xenopus Tropicalis* E3 Ubiquitin-Protein Ligase March2(41335) Protein, His-Tagged

Cat. No. RFL18276XF **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Recombinant Full Length <i>Xenopus tropicalis</i> E3 ubiquitin-protein ligase MARCH2(41335) Protein (Q28EX7) (1-246aa), fused to N-terminal His tag, was expressed in <i>E. coli</i> .
Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Source	<i>E. coli</i>
ProteinLength	Full Length (1-246)
Form	Lyophilized powder
AA Sequence	MTTGDCCHLPGLCDCTGSATFLKSLEESDLGRPQYVTQVTAKDGQLLSTVIKALGT QSD GPICRICHEGGNGERLLSPCDCTGTLGTVHKTCLEKWLSSSNTSYCELCHTEF AVERRPR PVTEWLKDPGPRNEKRTLFCDMVCFLFITPLAAISGWLCLRGAQDHLQF NSRLEAVGLIA LTIALFTIYVLWTLVSFRYHCQLYSEWRRTNQKVLILLIPDSKTAPTTH HSLSSKLLKSA SDETTV
Purity	Greater than 90% as determined by SDS-PAGE.
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Storage	Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

repeated freeze-thaw cycles.

Storage Buffer

Tris/PBS-based buffer, 6% Trehalose, pH 8.0

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

GENE INFORMATION

Gene Name

march2

Synonyms

marchf2; march2; TGas088c08.1; E3 ubiquitin-protein ligase MARCHF2; Membrane-associated RING finger protein 2; Membrane-associated RING-CH protein II; MARCH-II; RING-type E3 ubiquitin transferase MARCHF2

UniProt ID

[Q28EX7](#)

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA