

Recombinant Full Length Shigella Sonnei Nadh-Quinone Oxidoreductase Subunit K(Nuok) Protein, His-Tagged

Cat. No. RFL2061SF **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Recombinant Full Length Shigella sonnei NADH-quinone oxidoreductase subunit K(nuoK) Protein (Q3YZT1) (1-100aa), fused to N-terminal His tag, was expressed in E. coli.
Species	Shigella sonnei
Source	E.coli
ProteinLength	Full Length (1-100)
Form	Lyophilized powder
AA Sequence	MIPLQHGLILAAILFILGLTGLVIRRNLLFMLIGLEIMINASALAFVVGSYWGQTDGQV MYILAIISLAAAEASIGLALLQLHRRRQNLNIDSVSEMRG
Purity	Greater than 90% as determined by SDS-PAGE.
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Storage	Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.
Storage Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

GENE INFORMATION

Gene Name	nuoK
Synonyms	nuoK; SSON_2336; NADH-quinone oxidoreductase subunit K; NADH dehydrogenase I subunit K; NDH-1 subunit K
UniProt ID	Q3YZT1

 Tel: 1-631-559-9269 1-516-512-3133 Email: info@creative-biomart.com  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA