

Recombinant Full Length Chlamydia Muridarum 3-Deoxy-D-Manno-Octulosonic-Acid Transferase(Waaa) Protein, His-Tagged

Cat. No. RFL28215CF **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Recombinant Full Length Chlamydia muridarum 3-deoxy-D-manno-octulosonic-acid transferase(waaA) Protein (Q9PKI5) (1-430aa), fused to N-terminal His tag, was expressed in E. coli.
Species	Chlamydia muridarum
Source	E.coli
ProteinLength	Full Length (1-430)
Form	Lyophilized powder
AA Sequence	MRRWLTSRLYDAFLVAAFLAAAPRIFYKVVFHGKYINSWKIRFGVEKPQVKGEGPLV WFH GASVGEVSLLLEPLLKKWRQEFPDWRVVTACSEAGVYTAQRLYAPLGATV LPLDLSCI INPVVRSLSPPQVVFSEGDCWLHFLMGAKKLGAKAFLINGKLSSENSCKRF AFLKRLGRSY FAPDLLVLQDKVYKQRFMQIGIPEDKIQISGNLKTFIETETSINNRS WRKKLKLSPSD RLIVLGSMPKDVWADVAQHFNFSTKILWVPRHLEKLKEHAR LLEKAGISFGLWSKE DSLLQYDSLIVDAMGILKDLISAADLAFVGGTFDPLVGGHNLL EPLQKEVPLMFGPHIHS QSVLAELLRTKEVGVSVVDKENLLEAVENLLEDEKKRQAYI ERKGSFLKNAGTSFEHTWEI LKSQIACIKI
Purity	Greater than 90% as determined by SDS-PAGE.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Storage	Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.
Storage Buffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
GENE INFORMATION	
Gene Name	waaA
Synonyms	waaA; kdtA; TC_0480; 3-deoxy-D-manno-octulosonic acid transferase; Kdo transferase; Kdo(2-lipid IV(A 3-deoxy-D-manno-octulosonic acid transferase; Kdo-lipid IV(A 3-deoxy-D-manno-octulosonic acid transferase; Lipid IV(A 3-deoxy-D-manno-octulosonic acid tra
UniProt ID	Q9PKI5