

# Recombinant Full Length Campylobacter Jejuni Subsp. Jejuni Serotype O:2 Membrane Protein Insertase Yidc(Yidc) Protein, His-Tagged

**Cat. No.** RFL28534CF    **Lot. No.** (See product label)

## SPECIFICATION

|                         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|-------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b> | Recombinant Full Length Campylobacter jejuni subsp. jejuni serotype O:2 Membrane protein insertase YidC(yidC) Protein (Q9PNX7) (1-528aa), fused to N-terminal His tag, was expressed in E. coli.                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Species</b>          | Campylobacter jejuni subsp. jejuni serotype O:2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Source</b>           | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>ProteinLength</b>    | Full Length (1-528)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Form</b>             | Lyophilized powder                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>AA Sequence</b>      | <p>MNNSNNIFQQKRILLAVVISFLFFVIYDYFFIPKQPLKIEQNITQKNQQNTSINNTPNIQ<br/> NATTNTPSAALVSQDSVISKVQSKHFEAQIDSFGRISAFYLKDRKYQNEKGEFINLVS<br/> KE NSPYPLEMRFS DPSINSEAFKIPYVANASNLFVDENGSQLKLTQNL SGLKIEKDI<br/> TFYP KGNYEIEVKLSKNANYFISPGYRPNI AVDSYTVHGALVMDNKETIETYKDGDV<br/> EKDESAN NVVMTSAFD RYYATFFYNFDKPLNVAISKDANQNPIVFAYS DNEFKAGG<br/> YIGSKEHVILR SIDPRLEAVVEYGWFTFI AKPMFEFLNFLHQYIGNWGWAI VVMTLIV<br/> RIILFPLTYKSMI SMNKLKDLAPKMKDIRERYK GDPQKMNMHMMELYKKHGANPMS<br/> GCLPILLQIPIFFAIYR VLLNAIELKAAPWAFWIHDLSVMDPYFILPILMGATMFLQQLIT<br/> PMTIQDPMQAKIMKFL PVIFTFFFITFPAGLTLYW FVNNLCSLVQQWVINKIFAKEHH<br/> KKTSGA</p> |

 Tel: 1-631-559-9269    1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)     Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

|                         |                                                                                                                                                                                                                                                                                                                                                                                                             |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Purity</b>           | Greater than 90% as determined by SDS-PAGE.                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Notes</b>            | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Storage</b>          | Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.                                                                                                                                                                                                                                                                                              |
| <b>Storage Buffer</b>   | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Reconstitution</b>   | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>GENE INFORMATION</b> |                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Gene Name</b>        | yidC                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Synonyms</b>         | yidC; Cj0958c; Membrane protein insertase YidC; Foldase YidC; Membrane integrase YidC; Membrane protein YidC                                                                                                                                                                                                                                                                                                |
| <b>UniProt ID</b>       | <a href="#">Q9PNX7</a>                                                                                                                                                                                                                                                                                                                                                                                      |