

Recombinant *Pseudomonas taetrolens* ARGR protein, His-tagged

Cat. No. ARGR-285P **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Recombinant <i>Pseudomonas taetrolens</i> ARGR protein(I0J1I6)(393aa), fused with His Tag, was expressed in <i>E. coli</i> .
Species	<i>Pseudomonas taetrolens</i>
Source	<i>E. coli</i>
ProteinLength	393aa
Molecular Mass	42.4kDa
AA Sequence	MHHHHHHAPP LSMTDGVAQV NTQDSNAWVE INKAAFEHNI RTLQTALAGK SQICAVLKAD AYGHGIGLLM PSVIAMGVPC VGVASNEEAR VVRESGFKGQ LIRVRTAALS ELEAALPYNM EELVGNLDFA VKASLIAEDH GRPLVVHLGL NSSGMSRNGV DMTTAQGRRD AVAITKVPNL EVRAIMTHFA VEDAADVVRAG LKAFNQQAQW LMNVAQLDRS KITLHAANSF ATLEVPESHL DMVRPGGALF GDTVPSHTEY KRVMQFKSHV ASVNSYPKGN TVGYGRITYTL GRDSRLANIT VGYSDDGYRRA FTNKGIVLIN GHRVPVVGKV SMNTLMVDVT DAPDVKSGDE VVLFGHQGKA

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

EITQAEIEDI NGALLADLYT VWGNSNPIL KDQ

Endotoxin	< 1EU/μg
Purity	>95%(SDS-PAGE & RP-HPLC)
Storage	If the sample is used within 2-4 weeks, it can be stored at low temperature at 4°C. For long-term storage, it is recommended to add carrier protein (0.1% HSA or BSA) and store at -20~ -80°C. Avoid multiple freeze-thaws.
Reconstitution	It is recommended that sterile water be added to the vial to prepare a stock solution of 0.2 ug/ul. Centrifuge the vial at 4°C before opening to recover the entire contents.

 Tel: 1-631-559-9269 1-516-512-3133 Email: info@creative-biomart.com  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA