

ATP Agarose

Cat. No. ATP-001A **Lot. No.** (See product label)

SPECIFICATION

Product Overview	A ligand commonly used for this ATP Agarose resin is ATP (adenosine-5'-triphosphate) since a very large number of proteins such as kinases, motor proteins, and chaperones bind ATP with high affinity.
Cat.no	ATP-001A
Characteristics	ATP Agarose is used in the purification of ATP binding proteins. Any protein with an accessible ATP-binding domain will bind to it, regardless of activation state. The activated proteins will bind to the resin with higher affinity than the corresponding inactivated protein. Better recovery and yields may be obtained when the purification is followed by ion exchange chromatography.
Matrix	Highly cross-linked agarose
Average particle size	45 - 165 μ m
Ligand	ATP
Ligand concentration	2 – 5 ATP μ mol / ml drained gel
Chemical stability	Stable to all solutions commonly used in gel filtration including 8 M urea and 6 M guanidine hydrochloride. Not stable in organic solvents!
pH working range	7.5

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

pH CIP range

4 - 9

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA