

Recombinant Rat Mmp9 Protein Pre-coupled Magnetic Beads

Cat. No. Mmp9-3377R-B **Lot. No.** (See product label)

SPECIFICATION

Product Overview	The Recombinant protein was conjugated to magnetic beads. This ready-to-use, pre-coupled magnetic beads are in uniform particle size and narrow size distribution with large surface area, which is conducive to convenient and fast capture target molecules with high specificity and achieve magnetic separation. This product can be equipped with automation equipment for high-throughput operations.
Species	Rat
Source	HEK293
Form	Solution
Particle size	~2 µm
Beads Surface	Hydrophilic
Capacity	> 200 pmol rabbit IgG/ mg beads
Applications	Immunoassay, In vitro diagnostics, cell sorting, Immunoprecipitation/Co-precipitation, Protein/antibody separation and purification.
Stability	Stable for at least 6 months from the date of receipt of the product under proper storage and handling conditions.
Storage	2-8°C. Do not to freeze thaw the Beads

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Concentration	10mg beads/mL
----------------------	---------------

Storage Buffer	PBS buffer
-----------------------	------------

GENE INFORMATION

Gene Name	Mmp9 matrix metalloproteinase 9 [Rattus norvegicus]
------------------	---

Official Symbol	Mmp9
------------------------	------

Synonyms	MMP9; matrix metalloproteinase 9; matrix metalloproteinase-9; GELB; MMP-9; gelatinase B; 92 kDa gelatinase; 92 kDa type IV collagenase; 92-kDa type IV collagenase; matrix metalloproteinase 9 (gelatinase B 92-kDa type IV collagenase); matrix metalloproteinase 9 (gelatinase B, 92-kDa type IV collagenase);
-----------------	--

Gene ID	81687
----------------	-------

mRNA Refseq	NM_031055
--------------------	-----------

Protein Refseq	NP_112317
-----------------------	-----------

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA