

Recombinant Rhesus Macaque NFKBIL1 Protein Pre-coupled Magnetic Beads

Cat. No. NFKBIL1-2837R-B **Lot. No.** (See product label)

SPECIFICATION

Product Overview	The Recombinant protein was conjugated to magnetic beads. This ready-to-use, pre-coupled magnetic beads are in uniform particle size and narrow size distribution with large surface area, which is conducive to convenient and fast capture target molecules with high specificity and achieve magnetic separation. This product can be equipped with automation equipment for high-throughput operations.
Species	Rhesus macaque
Source	HEK293
Form	Solution
Particle size	~2 μ m
Beads Surface	Hydrophilic
Capacity	> 200 pmol rabbit IgG/ mg beads
Applications	Immunoassay, In vitro diagnostics, cell sorting, Immunoprecipitation/Co-precipitation, Protein/antibody separation and purification.
Stability	Stable for at least 6 months from the date of receipt of the product under proper storage and handling conditions.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Storage	2-8°C. Do not to freeze thaw the Beads
----------------	--

Concentration	10mg beads/mL
----------------------	---------------

Storage Buffer	PBS buffer
-----------------------	------------

GENE INFORMATION

Gene Name	NFKBIL1 nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor-like 1 [Macaca mulatta (Rhesus monkey)]
------------------	--

Official Symbol	NFKBIL1
------------------------	---------

Synonyms	NFKBIL1; NF-kappa-B inhibitor-like protein 1; I-kappa-B-like protein; ikappaBL; inhibitor of kappa B-like protein;
-----------------	--

Gene ID	715283
----------------	--------

mRNA Refseq	NM_001047144
--------------------	--------------

Protein Refseq	NP_001040609
-----------------------	--------------

UniProt ID	G7MRI1
-------------------	--------

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA