

Recombinant Human coronavirus HKU1 (isolate N5) (HCoV-HKU1) Spike Protein (S1+S2 ECD), His-tagged

Cat. No. Spike-547H **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Recombinant Human coronavirus HKU1 (isolate N5) (HCoV-HKU1) Spike Protein (S1+S2 ECD)(Q0ZME7.1)(Met1-Pro1295), fused with His tag, was expressed in Insect Cells.
Species	HCoV-HKU1
Source	Insect Cells
ProteinLength	Met1-Pro1295
Form	Lyophilized from sterile 20 mM PB, 300 mM NaCl, 10 % glycerol, pH 7.5Please contact us for any concerns or special requirements. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the hard copy of CoA.
Molecular Mass	The recombinant human coronavirus HKU1 (isolate N5) (HCoV-HKU1) spike protein (S1+S2 ECD, His Tag) consists of 1294 amino acids and predicts a molecular mass of 144.4 kDa.
Endotoxin	< 1.0 EU per µg protein as determined by the LAL method.
Purity	> 75 % as determined by SDS-PAGE.
Storage	Samples are stable for up to twelve months from date of receipt at -20°C to -80°C

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

It is recommended that sterile water be added to the vial to prepare a stock solution of 0.2 ug/ul. Centrifuge the vial at 4°C before opening to recover the entire contents.

GENE INFORMATION

Gene Name S surface glycoprotein [Severe acute respiratory syndrome coronavirus 2]

Official Symbol Spike

Synonyms Spike; S1 protein; Spike glycoprotein Subunit1; S glycoprotein Subunit1; Spike protein S1; COVID-19

Gene ID 43740568

Protein Refseq YP_009724390.1

UniProt ID P0DTC2

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA