

Active Spike S1 protein-coupled magnetic Beads

Cat. No. Spike-4921M **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Species	Sars-CoV-2
Source	HEK293
ProteinLength	Val 16 - Arg 685
Appearance	Powder mixture
Particle size	2 µm
Bio-activity	Immobilized 25 µg SARS-CoV-2 S1 protein to 1mg Beads, can bind the Anti-SARS-CoV-2 Spike S1 Antibody with an EC50 of 1.003 µg/mL.
Coupled amount of protein	≈293 pmol (23 µg) Spike S1/mg Beads
Capacity	>40 µg anti-SARS2-CoV-2 Spike S1 antibody/mg beads; >20 µg ACE2 protein/mg beads
Formulation	PBS, pH7.4, with 10% Trehalose
Reconstitution	2 mL sterile deionized water (1 mg beads/mL)
Background	The pre-coupled magnetic beads coupled with biotinylated SARS-CoV-2 Spike S1

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

protein to streptavidin conjugated magnetic beads, which can capture the Anti-SARS-CoV-2 antibody or ACE2 protein from cell or serum sample. The beads are in uniform size, narrow size distribution with large surface area and unique surface coating, which can help you get the best performance and highly reproducible results. This very first SARS-CoV-2 Spike protein S1-coupled magnetic beads will bring great convenience with minimum non-specific binding and developed protocols. This ready to use products could greatly save your time and hassle.

Application

This product is intended for immunocapture, biopanning and flow cytometry. This product is produced non-sterile.

Storage

Upon receipt, please store the lyophilized beads at -20°C for long term storage. The shelf life is 1 year at -20°C in lyophilized form.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA