

Active Spike RBD (B.1.351) protein-coupled magnetic Beads

Cat. No. Spike-4130M **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Species	Sars-CoV-2
Source	HEK293
ProteinLength	Arg 319 - Lys 537(K417N,E484K,N501Y)
Appearance	Powder mixture
Particle size	2 µm
Bio-activity	Immobilized 40 µg RBD protein /1mg beads can bind ACE2 with an EC50 of 2.395 µg/mL.
Coupled amount of protein	≈1445 pmol (40.75 µg) SARS-CoV-2 RBD protein/mg beads
Capacity	>300 pmol (33 µg) Human ACE2 ligand/mg beads
Formulation	PBS, pH7.4, with 10% Trehalose
Reconstitution	2 mL sterile deionized water (1 mg beads/mL)
Background	The SARS-CoV-2 Spike RBD (B.1.351) Coupled Magnetic Beads is produced by coupling biotinylated SARS-CoV-2 spike RBD to streptavidin-conjugated magnetic

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

beads. The spike RBD coupled to the beads contains the critical mutations (K417N, E484K, N501Y) identified in the South African variant (known as B.1.351 or 20C/501Y.V2). The pre-coupled beads are ready to use for capturing anti-SARS-CoV-2 antibody or ACE2 protein from your sample with high specificity.

Application

This product is intended for immunocapture, biopanning and flow cytometry. This product is produced non-sterile.

Storage

Upon receipt, please store the Beads at -20°C for 1 year in lyophilized state.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA