

GMP Genetically engineered DNA cytosine base editor

Cat. No. CBE-26G **Lot. No.** (See product label)

SPECIFICATION

Product Overview

CBE is a genetically engineered DNA cytosine base editor that efficiently modifies DNA within a range of 3-12 positions (with position No. 1 being farthest from the PAM site). By integrating the deaminase into the Cas enzyme, minimizes random binding to non-target DNA, reducing off-target effects, when combined with sgRNA (compatible with spCas9 sgRNA), forms an RNP complex. When the complex is transfected into cells and bound to the target site, it deaminates cytosine, creating uracil without causing DNA double-strand breaks. The cell's repair mechanism then changes the mutated site to thymine (T). This C->T mutation can knock out the target gene by generating a stop codon or mutating the alternative splicing site.

Source E.coli

Form 30 mM Tris, 300 mM NaCl, 0.1 mM EDTA, 1 mM DTT, 50% Glycerol, pH 8.0

Molecular Mass Greater than 190 KDa

Storage Store at -80 ±5°C. Avoid repeated freezing and thawing.

Reconstitution It is recommended that sterile water be added to the vial to prepare a stock solution of 0.2 ug/ul. Centrifuge the vial at 4°C before opening to recover the entire contents.

SEC-HPLC

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Flow cytometric analysis

