

## Recombinant MERS-CoV N protein, His-SUMO-tagged

**Cat. No.** N-3929V    **Lot. No.** (See product label)

### SPECIFICATION

|                         |                                                                                                                                                                               |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b> | Recombinant MERS-CoV N protein(AHI48804.1)(Met1-Asp413), fused with His and SUMO tag, was expressed in HEK293.                                                                |
| <b>Species</b>          | MERS-CoV                                                                                                                                                                      |
| <b>Source</b>           | HEK293                                                                                                                                                                        |
| <b>ProteinLength</b>    | Met1-Asp413                                                                                                                                                                   |
| <b>Form</b>             | Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.                                                                                                                   |
| <b>Molecular Mass</b>   | The protein has a calculated MW of 56.59 kDa. The protein migrates about 100 kDa under reducing condition (SDS-PAGE analysis).                                                |
| <b>Purity</b>           | >90% as determined by SDS-PAGE analysis.                                                                                                                                      |
| <b>Storage</b>          | Store at -20°C/-80°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles.                                                                |
| <b>Reconstitution</b>   | It is recommended that sterile water be added to the vial to prepare a stock solution of 0.2 ug/ul. Centrifuge the vial at 4°C before opening to recover the entire contents. |

 Tel: 1-631-559-9269    1-516-512-3133

 Email: [info@creative-biomart.com](mailto:info@creative-biomart.com)     Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA