

Recombinant TEV protease

Cat. No. Protease-532T **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Recombinant TEV protease was expressed in E. coli.
Species	TEV
Source	E.coli
Form	50mM PBNa, pH7.4, 50 mM NaCl, 5mM DTT, 1mM EDTA, 50% glycerol.
Molecular Mass	Approximately 51.5kD.
Endotoxin	<1EU/μg, determined by the LAL method.
Purity	≥95%.
Storage	rTEV in 50mM PBNa, pH7.4, 50 mM NaCl, 5mM DTT, 1mM EDTA, 50% glycerol should be stored at -20°C. It would be stable for 3 years when stored properly.
Reconstitution	It is recommended that sterile water be added to the vial to prepare a stock solution of 0.2 ug/ul. Centrifuge the vial at 4°C before opening to recover the entire contents.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA