

Blank Agarose (Control)

Cat. No. Agarose-002 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Blank Agarose
Sterilization	Not autoclavable!
Chemical stability	Stable to all solutions commonly used in gel filtration including 8 M urea and 6 M guanidine hydrochloride. Not stable in organic solvents!
pH stability	Short term: 4-9 Long term: 4-9
Maximum pressure	0.25 bar(3.6 psi)
Recommended linear flow rate	11.5 cm/h
Bead/Particle size	45-165 μ m
Applications	Suitable for control experiments with modified agarose.
Shelf Life	12 months after date of delivery
Storage	Store at 4 centigrade. Short term exposure (up to 1 week cumulative) to ambient temperature possible.
Storage Buffer	20 % Ethanol

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Shipping

Shipped at 4 centigrade.

Notes

For general laboratory use.

 Tel: 1-631-559-9269 1-516-512-3133 Email: info@creative-biomart.com  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA