

Gentamicin Detection Kit

Cat. No. KITF-048 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	The gentamicin is extracted from the fermentation medium of the genus Monosporinidae. It is an alkaline compound. It is a commonly used amino glycoside antibiotic and is widely used in culture medium preparation. Gentamicin quantitative detection kit adopts indirect competitive ELISA method, and the absorption value is negatively correlated with the content of Gentamicin in the sample. This kit only takes about an hour to operate, with a linear range of 0.1ng/mL-10ng/mL.
Size	96 tests/kit
Storage	Store in the dark at 2-8°C. Note: Opened but unused kits should still be stored in a dark place at 2-8°C.
Kit Components	Microtiter Plate; 96T(8 x 12 strips); 2-8°C Standard(100ng/mL); 1mL×1 tube; 2-8°C Sample Dilution; 30mL×1vial; 2-8°C Enzyme-labeled Secondary; 12mL×1vial; 2-8°C Antibody Antibody Working Solution; 7mL×1vial; 2-8°C Wash Buffer (20×); 30mL×1vial; 2-8°C Developing solution A; 8mL×1vial; 2-8°C(keep in dark) Developing solution B; 8mL×1vial; 2-8°C(keep in dark) Stop Solution; 15mL×1vial; 2-8°C

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA