

Biotinylated Bovine Serum Albumin

Cat. No. BSA-2939B **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Biotinylated Bovine Serum Albumin (Biotin-BSA) is prepared by covalent coupling of bovine serum albumin and Biotin. The combination is firm and difficult to be separated. It is generally used in combination with streptavidin and can be used in scientific research and in vitro diagnosis. In different application scenarios, it has the functions of signal amplification, signal bridge, signal shielding, etc.
Species	Bovine
Form	20mM PBS, pH7.4
Bio-activity	3-4mol Biotin/mol BSA
Storage	-20±5°C can be stored for 3 years; After thawing, it can be stored at 2-8°C for 2 weeks to avoid repeated freezing and thawing
Concentration	10mg/mL

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA