

Native DENV Serotype 2 Antigen

Cat. No. DEN-149 **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Inactivation	Formalin and detergent
Formulation	20 mM Phosphate Buffer, pH 7.4 with 140 mM NaCl and containing less than 0.1% CHAPS
Applications	ELISA
Storage	-70 centigrade.
Notes	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.
Type	Native
ClassID 1	Infectious Disease

GENE INFORMATION

Synonyms	DENV-Serotype 2
-----------------	-----------------

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA