

Recombinant Dengue Virus Serotype 1 VLP

Cat. No. DEN-169 **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Dengue Virus Serotype 1 VLP (Virus-Like Particle) is a new, unique product that has been developed in response to the need for highly immunogenic, concentrated Dengue antigens for use in Dengue focussed disease research and vaccine development (including use as an immunogen).

Background

Virus-Like Particles (VLPs) are an emerging vaccine technology. VLPs consist of protein shells comprising outer proteins specific to the virus in question. VLPs are more representative of how viral antigens are presented in vivo and while they are highly immunogenic they are non-infectious as they lack the core genetic material of the virus. Another important advantage of VLPs is that they more effectively activate key aspects of the immune response to achieve potent immune stimulation and to provide immunological memory. VLP-based vaccines have also shown to provide effective protection and are in current use for several diseases and in development for many others.

With more than one-third of the world's population living in areas at risk of transmission, dengue infection is a leading cause of illness and death in the tropics and subtropics (ref WHO for Guidelines for Diagnosis, Treatment, Prevention & Control, 2009) As many as 100 million people are infected yearly. Dengue is caused by any one of four related viruses transmitted by mosquitoes. There are no vaccines available to prevent infection with dengue virus. Dengue virus Serotype 2, Serotype 3 and Serotype 4 VLPs are also now available.

Purity

>95% pure by SDS-PAGE

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Sequence Strain	Puerto Rico/US/BID-V853/1998
Formulation	10mM sodium phosphate.20mM sodium citrate,154mM sodium chloride,pH7.4.
Stability	At +4 centigrade: 1 week. At -80 centigrade: 1 year.
Freezing	Can be frozen, but avoid multiple freeze/thaw cycles.
Storage	Short Term Storage: +2 to +8 centigrade Long Term Storage: -80 centigrade Avoid excessive mixing or shocking to prevent aggregation. Long term storage above -80 centigrade may result in aggregate formation.
Concentration	0.56 mg/mL
Notes	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.
Type	Recombinant
ClassID 1	Infectious Disease
GENE INFORMATION	
Synonyms	DENV Serotype 1 VLP

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA