

DNMT1 CoIP Agarose

Cat. No. DNMT1-01A **Lot. No.** (See product label)

SPECIFICATION


Product Overview	Immunoprecipitation of DNMT1 proteins with anti-DNMT1 Nanobody conjugated to agarose beads.
	• Fast, reliable & efficient one-step immunoprecipitation
	• Ready-to-use
	• No heavy & light antibody chains
	• Stable under harsh washing conditions
Specificity/Target	• Suitable for downstream mass spec analysis
Binding capacity	2-3 ug of DNMT1 protein per 10 uL bead slurry
Conjugate	Agarose beads; ~90 um (cross-linked 4% agarose beads)
Elution buffer	2xSDS-sample buffer (Lamml)
Type	Nanobody
Storage Buffer	20% ethanol
Applications	IP, Co-IP

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Western Blot



Tel: 1-631-559-9269 1-516-512-3133

Email: info@creative-biomart.com Fax: 1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA