

PARP1 CoIP Agarose beads, kit

Cat. No. PARP1-02A **Lot. No.** (See product label)

SPECIFICATION


Product Overview	The PARP1-Trap Agarose Kit contains PARP1-Trap Agarose, lysis, wash, and elution buffers for efficient Immunoprecipitation of p53. • Effectively immunoprecipitates PARP1 (human) • Functional after IP: on-bead parylation assay • No heavy & light chains in your downstream application
Specificity/Target	PARP1, Epitope: Within DNA binding domain of PARP1 does not bind to PARP2, 3, and 9
Tested Reactivity	Human
Conjugate	Agarose beads; bead size: ~ 90 µm (cross-linked 4 % agarose beads)
Elution buffer	SDS sample buffer 0.2 M glycine pH 2.5
Type	Nanobody
Storage Buffer	20% ethanol
Applications	IP, CoIP

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Western Blot



Tel: 1-631-559-9269 1-516-512-3133

Email: info@creative-biomart.com Fax: 1-631-938-8127

45-1 Ramsey Road, Shirley, NY 11967, USA