

Complement C3 assay kit (Immunoturbidimetry)

Cat. No. KITZ-076 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Based on the principle of immunoturbidimetry, the complement C3 antibody reacts with the complement C3 antigen to trigger an antigen-antibody reaction, forming an immune complex. The change in turbidity is measured at a wavelength of 340 nm, and the degree of change is proportional to the content of complement C3 in the sample.
Size	96 Samples
Storage	4°C
Shipping	Ice pack
Kit Components	Reagent One: Liquid, 25mL × 1 bottle, stored at 4°C.
	Reagent Two: Liquid, 5mL × 1 bottle, stored at 4°C.
	Standard Tube: Liquid, 0.1mL × 1 vial, stored at 4°C with a concentration of 2.68g/L.
Materials Required but Not Supplied	Microplate reader, 96-well plate, water bath, adjustable pipette, centrifuge, distilled water.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA