

Acetylcholintransferase (ChAT) activity kit

Cat. No. KITZ-088 **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Choline acetyltransferase (ChAT, EC 2.3.1.6) is the enzyme responsible for the synthesis of acetylcholine (Ach), regulating the metabolism of Ach. Ach is the primary excitatory neurotransmitter that regulates the movement of esophageal smooth muscle. Under the action of ChAT, acetyl-CoA and choline react to produce acetylcholine and coenzyme A. The product reacts with a chromogenic agent to yield an absorption peak at 412 nm, which is used to calculate the activity of ChAT. The enzyme-catalyzed reaction: $\text{acetyl-CoA} + \text{choline} = \text{CoA} + \text{O-acetylcholine}$

Size 48 Samples

Storage -20°C

Shipping Ice pack

Kit Components

Extraction Solution: Liquid, 60mL × 1 bottle, stored at 4°C.

Reagent 1: Powder, mg × 1 vial, stored at -20°C.

Before use, shake or centrifuge to bring the powder to the bottom, then add 1.2mL of distilled water to dissolve and prepare for use, which can be aliquoted and frozen at -20°C.

Reagent 2: Powder, mg × 1 vial, stored at 4°C.

Before use, shake or centrifuge to bring the powder to the bottom, then add 2.3mL of distilled water to dissolve and prepare for use, stored at -20°C.

Reagent 3: Liquid, 12mL × 1 bottle, stored at 4°C.

Reagent 4: Liquid, 4mL × 1 bottle, stored at 4°C.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

If solidified, it can be warmed in a 25°C water bath until completely dissolved before use.

Standard: Powder, 1mg × 1 vial, stored at 4°C. This reagent is used if a new calibration curve needs to be prepared.

**Materials Required
but Not Supplied**

Microplate reader, benchtop centrifuge, adjustable pipette, 96-well plate, mortar, ice, and distilled water

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA