

Endotoxin Removal Agarose resin

Cat. No. Endotoxin Removal-001A **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Endotoxin Removal beads is a product designed to remove endotoxins from biological protein products (including peptides, antibodies, polysaccharides, etc.). It involves attaching the modified polymyxin B to 4% agarose microspheres to specifically remove endotoxins. This process can reduce the endotoxin content in the sample to 0.1 EU/ml, and the sample recovery rate is high.
Matrix	4% agarose microspheres
Ligand	Modified polymyxin B
particle size	45-145 um
Maximum flow velocity	0.1 MPa, 1 bar
pH	5-10
Storage Buffer	1×PBS containing 20% ethanol
Storage	Store at 2-8°C.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA