

Endotoxin Removal Agarose resin, High Performance

Cat. No. Endotoxin Removal-001AHP **Lot. No.** (See product label)

SPECIFICATION

Product Overview	It is a filler used for removing endotoxins, with the ligand being polymyxin B. It has a high resolution and a high gravitational flow rate.
Matrix	4% highly cross-linked agarose
Ligand	Polymyxin B, 2mg/ml
particle size	30-70 um
Maximum flow velocity	300cm/h, 0.3 MPa
pH	3-10
Storage Buffer	20% ethanol, stored at 2-8°C, transported at 2-30°C.
Storage	Store at 2-8°C.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA