

Human Anti-Parvovirus B19(plasma) Reference standard

Cat. No. Parvovirus-174H **Lot. No.** (See product label)

SPECIFICATION

Product Overview	This material will serve as the primary biological standard for antibodies to Parvovirus B19.
Species	Human
Usage	No attempt should be made to weigh out any portion of the freeze-dried material prior to reconstitution. The contents of each ampoule should be reconstituted in 1ml of distilled water. Once reconstituted, diluted or aliquoted, users should determine the stability of the material according to their own method of preparation, storage and use.
Notes	This preparation is not for administration to humans. The preparation contains material of human origin, and either the final product or the source materials, from which it is derived, have been tested and found negative for HBsAg, anti-HIV and HCV RNA. As with all materials of biological origin, this preparation should be regarded as potentially hazardous to health. It should be used and discarded according to your own laboratorys safety procedures. Such safety procedures should include the wearing of protective gloves and avoiding the generation of aerosols. Care should be exercised in opening ampoules or vials, to avoid cuts.
Stability	Reference materials are held at Creative Biomart within assured, temperature-controlled storage facilities and they should be stored on receipt as indicated on the label. They remain valid with the assigned potency and status until withdrawn or amended.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA



Storage

Store ampoules at -20° or below until use Please note: because of the inherent stability of lyophilized material, Creative Biomart may ship these materials at ambient temperature.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA