

CM polystyrene resin

Cat. No. carboxy methyl-006P **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Creative Biomart provides series of ion exchange chromatography media for high efficient separation and purification of a wide range of biomolecules such as proteins, oligonucleotides, and peptides. The ion exchange resins are composed of rigid, spherical, highly cross-linked polystyrene/divinylbenzene or polyacrylic beads.
Cat.no	carboxy methyl-006P
Matrix	cross-linked polystyrene
Functional group	carboxy methyl group: $-O-CH_2COO-$
Protein	carboxy methyl
Dynamic binding capacity	50-80 mg/mL gel
Recommended flow rate	50-300 cm/h
Viscosity	The resin should be stored at 4-35 °C in 20 % ethanol
Chemical stability	Good chemical stability allows the use of the prepacked columns with most commonly used buffers and organic solvents, like for example 0.5 M NaOH, 50% ethanol, 20% isopropanol, 6 M GuHCl and 8 M urea.
Physical stability	Negligible volume variation due to changes in pH or ionic strength

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

pH working range	2-13
-------------------------	------

pH CIP range	1-14
---------------------	------

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA