

HBsAg Agarose

Cat. No. HBsAg-001A **Lot. No.** (See product label)

SPECIFICATION

Product Overview	HBsAg Agarose is widely used for the purification of Surface antigen of hepatitis B vaccine.
Cat.no	HBsAg-001A
Size	25 mL, 100 mL, 1 L
Matrix	Highly cross-linked agarose, 6%
Average particle size	90 µm
Ligand	HBsAg
Ligand concentration	10-40 µmol / ml wet gel
Recommended flow rate	<150>
Viscosity	store at 4-30°C
Chemical stability	The quality of resin doesn't change at all when dip in 1 M NaOH, 0.01 M HCl, 6 M guanidine hydrochloride and 8 M urea for 7 days under room temperature.
pH working range	3-13

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

pH CIP range	2-14
---------------------	------

Temperature Stability	4-30 °C
----------------------------------	---------

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA