

Mouse Anti-Rabbit IgG beads

Cat. No. mouse anti-rabbit IgG-17A **Lot. No.** (See product label)

SPECIFICATION

Product Overview

This product contains sepharose beads coated with mouse anti-rabbit IgG. Mouse anti-rabbit IgG is affinity purified and conjugated to the beads at 1-1.5 mg/ml ratio. The beads may be used for purification of rabbit IgG class of antibodies or proteins. The beads may also be used for immunoprecipitation (IP). For a complete use guide, please follow the protocol section. When used and stored properly, the beads may be used repeatedly for up to 100 times.

Cat.no mouse anti-rabbit IgG-17A

Matrix Highly cross-linked agarose, 4%

Functional group mouse anti-rabbit IgG

Average particle size 45-165 μ m

Ligand Mouse Anti-Rabbit IgG beads

Ligand concentration 1-1.5 mg/ml

Dynamic binding capacity 4-7 mg IgG per ml

Recommended flow rate 1-3 ml/min/cm²

Temperature 4°C to 25°C

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Stability

Elution buffer	2% sucrose, adjusted pH to 2-3 by CH ₃ COOH
-----------------------	--

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA