

TargetLock Protein G Beads

Cat. No. Protein G-34A **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Recombinant Protein-G used in our Protein G beads is a genetically engineered form of Protein-G. Non-essential regions have been removed while leaving the IgG binding sites intact. The recombinant Protein G has been immobilized on Sepharose beads, using our Biotech's TargetLock chemistry, in order to create a versatile product for both the purification of various mammalian immunoglobulins and for isolating immune complexes by immunoprecipitation (IP).
-------------------------	---

Cat.no	Protein G-34A
---------------	---------------

Matrix	Highly cross-linked spherical agarose, 4%
---------------	---

Average particle size	40-165 μ m
------------------------------	----------------

Dynamic binding capacity	~4-8 mg rabbit IgG per ml of beads
---------------------------------	------------------------------------

Recommended flow rate	1-3 ml/min/ cm ²
------------------------------	-----------------------------

pH working range	2-11
-------------------------	------

Elution buffer	0.1M glycine buffer Adjust pH to 2.8 by 1N HCl
-----------------------	--

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA