

HTP resin

Cat. No. Hydroxyapatite-49A **Lot. No.** (See product label)

SPECIFICATION

Product Overview	The resin is dried by a unique process developed by us Creative Biomart. It may be stored without refrigeration. When suspended in buffer, it has characteristics similar to Creative Biomart HT resin; 1 g hydrates to 2–3 ml.
-------------------------	---

Cat.no	Hydroxyapatite-49A
---------------	--------------------

Matrix	Highly cross-linked agarose
---------------	-----------------------------

Recommended flow rate	35-100 cm/h
------------------------------	-------------

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA