

Annexin V-FITC/7-AAD Apoptosis Kit

Cat. No. Kit-0096 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Annexin V-FITC/7-AAD Apoptosis Kit is used for apoptosis detection.
Description	Annexin A5 (or Annexin V) is a cellular protein in the annexin group. Annexin A5 is used as a probe to detect cells that have expressed phosphatidylserine (PS) on the cell surface, an event found in apoptosis as well as other forms of cell death. The annexin A5 affinity assay typically uses a conjugate of annexin V and a fluorescent or enzymatic label, biotin or other tags, or a radioelement, in a suitable buffer (annexin V binding to PS is Ca²⁺ dependent). The assay combines annexin V staining of PS membrane events with the staining of cell nucleus with PI or AAD-7 in living cells to distinguish from dead cells. Annexin V Apoptosis Detection Kit is based on the observation that soon after initiating apoptosis, cells translocate the membrane phosphatidylserine (PS) from the inner face of the plasma membrane to the cell surface. Once on the cell surface, PS can be easily detected by staining with a fluorescent conjugate of Annexin V, a protein that has a high affinity for PS. Detection can be analyzed by flow cytometry or by fluorescence microscopy.
Usage	For research use only (RUO)
Storage	Store the kit at 4°C.
Kit Components	Annexin V-FITC: Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze. 500 µL 7-AAD Viability Staining solution: Store undiluted at 4°C and protected from prolonged exposure to light. Do not freeze. 500 µL Binding buffer (10X): Store undiluted at 4°C. For long term storage, aliquot and store at -20°C. 30

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

	mL
Detection method	Fluorimetric
Compatible Sample Types	Cell

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA