

Cell Viability Green Fluorescence Assay Kit

Cat. No. Kit-0200 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Cell Viability Assay Kit (Green Fluorescence) is a fluorescent assay used to label cells for fluorescence microplate and microscopic investigations of cellular functions.
Description	There are a variety of parameters that can be used to monitor cell viability. The proprietary green fluorescent dye used in the kit is a hydrophobic compound. It easily permeates intact live cells and gets enhanced fluorescence upon entering into live cells. The hydrolysis of the non-fluorescent substrate by intracellular esterases generates a strongly green fluorescent hydrophilic product that is well-retained in the cell cytoplasm. The green fluorophore generated by the non-fluorescent substrate used in the kit has the spectral properties of fluorescein at Ex/Em = ~490 nm/520 nm. When well excited with the Argon Laser at 488 nm, the fluoreophore emits intense green fluorescence at ~520 nm.
Applications	The kit provides all the essential components with an optimized cell-labeling protocol for fluorescence microplate assays.
Usage	For research use only (RUO)
Storage	Store the kit at -20°C and avoid from light.
Kit Components	Component A: CytoCalcein Green (lyophilized) 5 vials Component B: DMSO 200 µL Component C: Assay Buffer 50 mL
Detection method	Fluorometric

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

**Compatible Sample
Types**

Cells

Features & Benefits

Robust: Higher maximum signal with lower variation across the plate. Convenient:
Formulated to have minimal hands-on time. Rapid Dye Loading: Dye loading at RT for
30 min to 1 hr. Versatile Applications: Compatible with many cell lines and targets.

 Tel: 1-631-559-9269 1-516-512-3133 Email: info@creative-biomart.com  Fax: 1-631-938-8127 45-1 Ramsey Road, Shirley, NY 11967, USA