

Choline/Acetylcholine Quantification Kit

Cat. No. Kit-0226 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Choline/Acetylcholine Quantification Kit is used for quantifying choline and acetylcholine by either a colorimetric or fluorometric method.
Description	Choline and acetylcholine play important roles in many biological processes. The Choline/Acetylcholine Quantification Kit provides a simple and sensitive means for quantifying Choline and Acetylcholine by either a colorimetric or fluorometric method. In the assay free choline is oxidized to betaine, via the intermediate betaine aldehyde. The reaction generates products which react with the Choline Probe to generate color ($\lambda = 570$ nm), and fluorescence (Ex/Em 535/587 nm). Acetylcholine can be converted to choline by adding acetylcholinesterase to the reaction. The kit can detect choline and acetylcholine (total choline – free choline) in various biological samples such as in blood, cells, culture media, fermentation media, etc. There is no need for pretreatment or purification of samples. The kit can detect 10 pmol~5 nmol of choline or acetylcholine.
Target Species	Mammals
Usage	For research use only (RUO)
Storage	Store kit at -20°C. Warm the assay buffer to room temperature before use. Briefly centrifuge vials before opening. Read the entire protocol before performing the assay.
Kit Components	Choline Assay Buffer 25 mlCholine Probe* 200 μ lCholine Enzyme Mix 1 VialAcetylcholinesterase 1 VialCholine Standard (5 μ mol) 1 VialNote: *The fluorometric

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

assay is ~10 times more sensitive than the colorimetric assay. Use 0.4 μ l of the probe per reaction to decrease background/increase detection sensitivity significantly.

Detection method	Colorimetric, Fluorometric
Compatible Sample Types	Biological Fluid, Cell Culture Supernatant, Plasma, Serum, Tissue Culture Supernatant, Urine

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA