

Cyclooxygenase 1 Inhibitor Screening Fluorometric Kit

Cat. No. Kit-0256 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Cyclooxygenase 1 (COX1) Inhibitor Screening Kit (Fluorometric) offers a rapid, simple, sensitive, and reliable test suitable for high-throughput screening of COX-1 inhibitors. The assay is based on the fluorometric detection of Prostaglandin G2, the intermediate product generated by the COX enzyme.
Size	100 tests
Description	Cyclooxygenase (COX), also known as prostaglandin-endoperoxide synthase (PTGS, EC 1.14.99.1), is an enzyme that is responsible for the formation of important biological mediators called prostanoids, including prostaglandins, prostacyclin and thromboxane. COX is the central enzyme in the biosynthetic pathway to prostanoids from arachidonic acid. There are two known isoenzymes: COX-1 and COX-2. COX-1 is constitutively expressed in many tissues and is the predominant form in gastric mucosa and in kidney. COX-2 is not expressed under normal conditions in most cells, but elevated levels are found during inflammation. Pharmacological inhibition of COX by nonsteroidal anti-inflammatory drugs (NSAID) can provide relief from the symptoms of inflammation and pain.
Applications	Inhibition Assaymore details
Target Species	Reacts with: Species independent
Storage	Store at -20°C.
Kit Components	Components: 100 tests; Arachidonic Acid: 1 x 50µl COX Assay Buffer: 1 x 25ml; COX

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Cofactor (in DMSO): 1 x 20 μ l COX Probe (in DMSO): 1 x 200 μ l COX-1, Ovine: 1 vial;
NaOH: 1 x 500 μ l SC560, COX-1 inhibitor (in DMSO): 1 x 100 μ l

Detection method Fluorescent

Compatible Sample Types Inhibitor compounds

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA