

Galactose and Lactose Assay Kit

Cat. No. Kit-0356 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Galactose and Lactose Assay Kit is used for measuring of galactose and lactose levels.
Description	Galactose (C ₆ H ₁₂ O ₆ FW: 180.16) and lactose (C ₁₂ H ₂₂ O ₁₁ FW: 342.3) are naturally occurring sugars. Lactose consists of one galactose and one glucose. Some people, particularly infants, lack the enzyme necessary to digest galactose leading to galactose accumulation in blood (Galactosemia) causing enlarged liver, renal failure, cataracts and brain damage. Galactose and Lactose Assay Kit provides a simple, convenient means for direct measurement of galactose and lactose levels in various biological samples (serum, plasma, other body fluids, food, growth media, etc.). To detect Galactose, galactose oxidase specifically oxidizes free galactose, generating a product that reacts with the Galactose Probe to produce color (O.D 570 nm) and fluorescence (Ex/Em 535/587). To detect lactose, lactose is hydrolyzed by lactase to generate free galactose which is then measured. Thus, Lactose level = Total Galactose - Free Galactose.
Target Species	Mammals
Usage	For research use only (RUO)
Storage	Store kit at -20°C, protect from light. Allow reagents warm to room temperature and briefly centrifuge vials before opening.
Kit Components	Galactose Assay Buffer 25 mLGalactose Probe* (DMSO solution) 0.2 mLLactase

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

(Lyophilized) 1 vialGalactose Enzyme Mix (Lyophilized) 1 vialHRP (Lyophilized) 1 vialGalactose Standard (100 nmol/μL) 100 μL

Detection method Colorimetric, Fluorometric

Compatible Sample Types Biological Fluid, Cell Culture Supernatant, Plasma, Serum, Tissue Culture Supernatant, Urine

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA