

HDAC Activity Assay Kit

Cat. No. Kit-0422 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	The HDAC Activity Assay Kit is designed to measure HDAC activity in cell or nuclear extracts, immunoprecipitates or purified enzymes.
Size	1 kit/96 Tests
Usage	The Calbiochem HDAC Activity assay is a complete assay system designed to measure histone deacetylase (HDAC) activity in cell or nuclear extracts, immunoprecipitates, or purified enzymes. It has been used successfully with preparations of all class I HDACs-HDAC1, HDAC2, HDAC3 and HDAC8-with class II HDACs 4-7, 9 and 10 and with yeast Sir2 and its human homolog, SIRT1.
Storage	Storage:< -70°C Storage Conditions: Upon arrival store the HeLa Cell Nuclear Extract, HDAC Substrate, HDAC Developer Concentrate, Trichostatin A (HDAC Inhibitor), Deacetylated Standard, and HDAC Assay Buffer at -70°C. Store both of the 1/2 Volume Microplates at room temperature.Avoid freeze/thaw: Avoid freeze/thaw Do not freeze: No
Compatible Sample Types	cell extracts, immunoprecipitates, or purified enzymes

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA