

JC-1 Mitochondrial Membrane Potential Assay Kit

Cat. No. Kit-0487 Lot. No. (See product label)

SPECIFICATION

Product Overview	JC-1 Mitochondrial Membrane Potential Assay Kit is a fluorescence-based method for tracing the behavior of mitochondria.
Description	$\Delta\Psi$ is an important parameter of mitochondrial function and has been used as an indicator of cell health. Variation of $\Delta\Psi$ have been previously studied using cationic dyes such as rhodamine-123 (Rh123) and DiOC6. More recently, a cytofluorimetric, lipophilic cationic dye, 5,5',6,6'-tetrachloro-1,1',3,3'-tetraethylbenzimidazole carbocyanine iodide (JC-1), has been developed. JC-1 has advantages over other cationic dyes in that it enters the mitochondria and changes its fluorescent properties based on the aggregation of the probe. In healthy cells with high $\Delta\Psi$, JC-1 forms complexes known as J-aggregates with intense red fluorescence. However, in cells with low $\Delta\Psi$, JC-1 remains in the monomeric form, which exhibits green fluorescence. The higher the ratio of red to green fluorescence, the higher the polarization of the mitochondrial membrane.
Usage	For research use only (RUO)
Storage	JC-1 Reagent -20°C Cell-Based Assay Buffer Tablet Room Temperature
Kit Components	JC-1 Reagent 500 μ L Cell-Based Assay Buffer Tablet 3 tablets

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA