

PLTP Inhibitor Drug Screening Kit

Cat. No. Kit-0710 **Lot. No.** (See product label)

SPECIFICATION

Product Overview	Simple procedure; takes ~40 minutesFast and convenientProvided all necessary buffers and reagents, as well as positive control for accurate assay of activity of plasma phospholipid transfer protein (PLTP).
Size	100 assays
Description	Plasma phospholipid transfer protein (PLTP) is thought to play a major role in the facilitated transfer of phospholipids between lipoproteins and in the modulation of high-density lipoprotein (HDL) particle size and composition. PLTP-facilitated lipid transfer activity is related to HDL and LDL metabolism, as well as lipoprotein lipase activity, adiposity, and insulin resistance. The PLTP Drug Screening Kit utilizes rabbit serum as PLTP activity source to screen PLTP inhibitors. The assay uses a donor molecule containing a fluorescent self-quenched phospholipid that is transferred to an acceptor molecule in the presence of PLTP. PLTP-mediated transfer of the fluorescent phospholipid to the acceptor molecule results in an increase in fluorescence (Ex/Em: 465/535). Inhibitor of PLTP will inhibit the lipid transfer and therefore decrease fluorescence intensity. The kit provides sufficient reagents for 100 PLTP inhibitor screening assays.
Target Species	Mammalian
Usage	This product is furnished for LABORATORY RESEARCH USE ONLY. Not for diagnostic or therapeutic use.

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

Kit Components	PLTP Donor MoleculePLTP Acceptor MoleculePLTP Assay Buffer (10X)Positive Control (Rabbit Serum)
Detection method	Fluorometer
Compatible Sample Types	Cell Culture SupernatantsPlasmaSerumUrine

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA