

Thioredoxin-1 assay Kit

Cat. No. Kit-0832 **Lot. No.** (See product label)

SPECIFICATION

Product Overview

Size 2x96 wells

Description

The Thioredoxin-1 assay kit provides a simple, homogenous assay for screening modulators of Trx-1 enzymatic activity in microplates. This is accomplished by monitoring the Trx-1 catalyzed reduction of insulin in the presence of Dithiothreitol (DTT), resulting in the formation of insulin aggregates which then bind avidly to the red-emitting fluorogenic ProteoStat Trx-1 detection dye. The fluorescence-based assay provides a vastly improved assay signal window, improved lower detection limit, and superior Z'-score (>0.8) relative to turbidity measurements. Intra-plate and inter-plate CVs using this assay are typically 3-6%. Human recombinant thioredoxin-1 enzyme and the thioredoxin inhibitor, aurothiomalate, are included in the kit, as well as all reagents necessary for monitoring changes in Trx-1 activity. The ProteoStat Thioredoxin-1 assay kit is capable of providing a quantitative readout of Trx-1 enzymatic activity in a robust and high-throughput fashion and can be applied to identification of Trx-1 inhibitors from chemical libraries.

Applications

This kit provides a simple, homogenous assay for screening modulators of Thioredoxin-1 enzymatic activity.

Storage

-80°C

Kit Components

ProteoStat Trx-1 detection reagent: 20µL
Thioredoxin-1 (human) (recombinant): 2 vials x 80µL
Insulin: 2 vials x 1.8µmoles
Aurothiomalate (inhibitor): 2 vials x

 Tel: 1-631-559-9269 1-516-512-3133

 Email: info@creative-biomart.com  Fax: 1-631-938-8127

 45-1 Ramsey Road, Shirley, NY 11967, USA

25nmolSPBE buffer: 25mlStop reagent: 1mlDTT: 2 vials x 1.3mlDeionized water: 5ml